R.R.Gen

Model 1.6.1

(Raymond Rife Generator)

Access Point:	R.R.Gen
	192.168.4.1
Password:	12345678

User Manual

Premise

Everything created has its own vibration.
According to quantum mechanics, matter is never inert, but is constantly in a state of motion, of continuous vibration.
The Austrian physicist Fritjof Capra said:

«Each particle perpetually sings its song».

Everything that makes up reality vibrates.
Many scientists have studied the healing effect of frequencies.
Dr. Royal Raymond Rife in the 1920s developed a device based on frequencies to heal almost any disease.
Even inanimate and dense objects such as stones that appear to us as solid matter are, in fact, forms of energy that vibrate, albeit at very slow frequencies.

Everything in the Universe is vibrating energy and therefore sound.

We are no longer used to listening, to perceiving these vibrations.
The main difference between our civilization and some of the advanced ancient civilizations is precisely in the concept of energy.
Ours is based on the production of energy from various sources: coal, oil and gas (dirty, polluting and expensive), nuclear (contaminating and extremely expensive).
The concept of energy of the ancients was instead based on the amplification of the energy produced by natural sources through the use of geometric shapes, such as the pyramid, the cone and the circle, and the elements of sacred geometry such as the equilateral triangle, the hexagon, the π , the ϕ or Golden Section and the Flower of Life.

R.R.Gen

Model 1.6.1

(Raymond Rife Generator)

All these geometric shapes and those generated by the Golden Section ratio cause the movement of energy and increase its intensity.

They amplify the natural energy sources existing underground produced by metals and minerals, by underground waterways, by quartz crystals.

Every chemical reaction, every element, atom or molecule is held together by small electromagnetic bonds created by resonance and frequency.

There is a basic symphony, a primordial vibration that guides everything that is present in the Universe and every living being resonates with this symphony.

Reprogramming cells with frequencies

This technology has numerous fields of application to improve our existence on our beautiful planet.

Operation

The R.R.Gen generator is a device that generates pulsed electromagnetic impulses (PEMF) that help balance the energy and life force of the human body, for greater overall well-being.

With normal use, you can experience good results and general well-being.

The power generated is minimal and has no side effects, even if the wrong treatment is chosen.

User Manual

- **Turning it on and off**
 - To turn the device on or off, operate the slide switch referring to the silkscreen on the panel
- **Connecting to your smartphone**
 - Turning on the device
 - Open the WiFi management menu


R.R.Gen

Model 1.6.1

(Raymond Rife Generator)

- Turn on the WiFi and search for the one named "R.R.Gen" followed by a number (this is a unique number that is different for each device) and select it
- in most phones the main page should open automatically, otherwise you will need to pull down to refresh or tap on the notification that opens after connecting
- if you still have difficulty connecting, Perform the procedure described in the next point, which is required for use with a PC
- For use with a PC (with WiFi connection) open a browser and search for the page at 192.168.4.1 and the main menu page will open
- **Navigating the menus**
 - The main page or MENU has buttons for navigation:
 - selection of treatment
 - help / quick guide
 - a blue button for settings
- After a selection or execution of a treatment, at the beginning, there may be green buttons to link to "favorite" treatments (maximum 3)
- **The settings**
 - at the top is displayed:
 - the version and date of the last software update
 - the name of the active treatment package
 - the serial number of the device (3 groups of 6 digits) that is required for the purchase of other packages with additional treatments in addition to those already present
 - a series of possible choices
 - the name of the user (it is possible to enter up to 5 different users that can be used to define groups of personalized treatments)
 - selections for preferred language
 - selection for active treatment package
 - the number of repetitions of the treatments to be performed
 - the intensity of the green LED light (for night or in bright light)
 - the cancel and confirm buttons

- **Searching for and selecting a treatment**
 - after selecting a treatment, you can perform several actions:
 - perform ONLY the selected treatment once (red button)
 - if you have previously entered a number of repetitions, perform the selected treatment and possibly the others in your favorites
 - simply add the selected treatment that will be placed in the first position
- **The execution of a treatment**
 - you can proceed in 2 ways:
 - you can perform a treatment starting from the main page, choosing it from the favorites and confirming its execution
 - you can search for it among those present in the active package, select it and confirm its execution
 - in any case the chosen treatment will be added to the first position of the favorites
 - the duration of the treatments depends on the number of frequencies of the treatment itself and the execution times for each frequency (on average each treatment lasts approximately 30 minutes)

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